

Novel Approach to RCC Early Diagnosis and Therapeutic Monitoring using Volatile Organic Compounds

VIVEK NARAYAN, MD MSCE

UNIVERSITY OF PENNSYLVANIA

GU MEDICAL ONCOLOGY

Clinical and Translational Research for Renal and Prostate Cancers

Early Phase Clinical Trials / Novel Therapeutics

Investigator Initiated Trials

- Phase 1b Avelumab / Gemcitabine for sarcomatoid mRCC
- Phase II Maintenance Niraparib for Biomarker Selected mCRPC
- Phase 1 First-in-Human PSMA Re-directed / TGFb-resistant CART for mCRPC
- Randomized neoadjuvant IO v. IO/TKI in high-risk RCC

Industry & Cooperative Group Trials

Cardiac Toxicity in GU Malignancies

- Multi-Institutional Prospective Study of Sunitinib-Induced Cardiac Toxicity
- Pilot Cohort Study of Myocardial Blood Flow Reserve using Rb82 Myocardial PET in Prostate Cancer
- Two Multi-Institutional Randomized Trials of CV Risk Mitigation Strategies
 - Intensive BP control vs Standard Care in mRCC
 - Modified Cardiac Rehab program in PCa

Translational / Biomarker

- Peripheral Blood Predictive Biomarkers in mRCC receiving anti-PD1 therapy
- Optimization / Validation and Early Clinical Testing of Novel Diagnostic/Classifier for Localized RCC

Brief Background:

- Up to 25% of incident RCC presents with locally advanced or metastatic disease
 - Initial clinical staging significantly impacts long-term clinical outcomes
- Critical need for early diagnostic, as well as improved methods for post-nephrectomy surveillance and therapeutic monitoring in advanced disease
- No current tools exist to efficiently conduct targeted screening in an average-risk population


Brief Background:

- Concept inspired by natural abilities of mammalian olfaction
 - Capitalizes on differential production of Volatile Organic Compounds (VOCs / 'odorants') by tumor cell metabolism
 - Controlled studies demonstrate sensitivity of canine detection in variety of localized malignancies
- Potential untapped target for early detection, surveillance, therapeutic monitoring?






CONCEPT OVERVIEW of VOC Detection

- Hypothesis: VOCs from healthy cells differ from benign tumor cells which differ from malignant cells
- Sensors have overlapping selectivity and are not specific to a chemical compound
- Electronic signals are analyzed using pattern recognition via machine learning to build a classification algorithm



Assay Design:

- Carbon nanotubes (CNT) connect the source and drain electrodes of a single field effect transistor (FET)
- CNTs are functionalized with single stranded DNA which increases response to VOCs by introducing binding sites
- Each CNT-FET is functionalized with the same DNA oligomer to form a unique odorant sensor type
- Sensor array comprises 10 CNT-FETs, yielding 10 sensor types per array
- VOC analytes interact with the DNA, altering the electrical properties of the FET, resulting in measurable shift in the current




Preliminary Experience with Ovarian Cancer:

Training Set: 24 plasma samples.

~ 8 each: Malignant (1 early stage), Benign, Control

Validation Set: 34 plasma samples.

11 Malignant (**5 early stage**), 11 Benign, 12 Control



Using Machine Learning for Classification and Statistical Validation:

(N=58 samples)

Stratified K-Fold Cross-validation

Success Rate: 95%

Early Stage Rate: 100%

False Positive: 0%

False Negative: 5%

		CLASSIFIED		
		Cancer	Benign	Healthy
ACTUAL	Cancer	10	1	0
	Benign	0	8	3
	Healthy	0	0	12

Preliminary Aims:

AIM 1: Evaluate for detection of VOC signature using DNA-NT sensor array technique in RCC patient samples.

- Further refine the DNA-NT sensor array for RCC by screening against large number of patient samples (from benign, healthy control, localized and metastatic disease settings).

AIM 2: Validate and benchmark the VOC sensor array

- Assessment of reproducibility, classification accuracy, population differences

AIM 3: Preliminary assessment of assay in clinical testing

- Presurgical clinical staging – correlation with T-stage, histopathology, nuclear grade
- Potential utility for post-nephrectomy surveillance